

# **Aflatoxion B1 Test Kit**

Cat No.:T19016

#### Introduction

It is based on the colloidal gold competition method to detect the residues of Aflatoxion B1 in food, feed ingredients and some finished feeds through specific antigen-antibody reaction. It is suitable for various enterprises, testing institutions and supervision departments detection.

The Detection Limit: 5~10ppb

## Can detect sample

Grain, Feeding, Parts of finished Feeding

## **Kit Components**

- 1. Test Strips, 12 pack, each pack 8 strips, 8pcs red reagent
- 2. 8pcs white microwells/pack
- 3. Manual Instruction

# Components Required, But Not Provided

- 1. Weight Scale (Accurancy above 0.1g)
- 2. 200µL Pipette
- 3. 50ml centrifuge tube
- 4. Ethanol (analytical grade), pure water
- 5. Vortex instrument, Timer

## Storage

At 4~30°C, out of sun light & moisture, sealed, do not freeze.

## **Solution Prepare**

50% ethanol-water solution: 50mL ethanol (analytical grade) + 50mL pure water

## **Sample Preparation**

Prepare the sample, and the sample should be in room temperature, no need to heat or refrigerate.


- 1. Weight 0.5g±0.05g grain sample into centrifuge tube. (Strong water absorption feed, such as Bran, Cotton aphid, Cottonseed etc. Weight 2.5g sample)
- 2. Add sample into centrifuge tube, mixing with 20ml 50% methanol Pure water.
- 3. Shake up and down the tube for 2 mins
- 4. Wait for 10mins, take the supernatant liquid. The Supernatant liquid is the sample liquid.

#### **Test Procedure**

- Take out the strips for assay. The unused test kit must be sealed and out of sunlight and moisture. Let the sample liquid turn to 20°C-25°C.
- 2. Drop 200ul sample liquid into the white microwell, Use pipette pietting up and down at least 6times, mix the liquid completely.










- 3. Then transfer 200ul mix sample liquid from the white microwell into the red microwell, pipetting up and down at least 6times, mix the liquid completely.
- 4. Timer 5mins at room temperature.
- 5. Insert the strip into microwell, with BALLYA Logo handle end is up.
- 6. Wait for 5mins, take out the sample pad. Reader the results.

# **Interpretation of Test Results**

A. Reading the Results

Results must be read within 5 minutes.



Negative: T line stronger than C line.

Positive: Other results are positive.

Invalid: Only T line appear or no C line.

**Performance Information** 

LOD: 1mg/kg (mg/kg=ppm)

#### **Precautions and Warnings**

- 1. Make sure let the strips turn to room temperature then for use.
- 2. Use the strips in the shelf time.
- 3. Keep the kits in sealed when do not use.
- 4. The kit is single-use, do not reuse.
- 5. This method is a screening method. If there is a positive sample, please confirm it by confirmatory method.

#### **Shelf Life**

12 months

# **Ballya International Limited**

Add: 5th Floor, No.241, Shangye Rd, Huadu, Guangzhou, China. 510810

Tel: 86 20 39479163 <u>www.ballyabio.com</u> info@ballyabio.com